


RETRACTION REVOLUTION

INTRODUCING THE **VIEWSITE™**,
THE NEXT GENERATION OF BRAIN RETRACTION.



The **ViewSite™** Brain Access System provides excellent surgical access to the operative site, allows for superior binocular vision of the surgical area, distributes the brain evenly upon insertion allowing the surgeon a seamless entry to the target, and eliminates the need to pull the brain tissue in any one direction. Its unique design allows for maximal lighting without the need for any fiberoptics. It is compatible with



most surgical arms and may be used with neuronavigational systems.

The **ViewSite™** Brain Access System is available in various shapes and sizes to allow for a diverse surgical application. **The ViewSite™ TC Model is available in** 3cm, 5cm and 7cm lengths and 12 mm, 17mm, 21mm and 28 mm widths. **The EC Model line includes** 3cm, 5cm and 7cm lengths and 34mm widths.

www.vycormedical.com



Targeting Solutions in Neurosurgery.